


Charles B.

SCHWED

Co., LLC

CATALOG



**ENDODONTIC
INSTRUMENTS AND
OTHER SPECIALIZED
DENTAL PRODUCTS**

1 Padanaram Rd, Ste 110, Peacock Alley, Danbury, CT, 06811

Tel: 475-289-3197

www.schwed.com

**WELCOME
TO THE
FAMILY,
WHERE
YOU
ARE
#1!**

Table of contents

Terms (2)

Stainless Steel Instruments (3-15)

K - Files (4-5)

Hedstrom Files (6-7)

Reamers (8-9)

C - Files (10)

Flexicut Files (11)

Barbed Broaches (12)

Finger Pluggers (13)

Finger Spreaders (13)

Gates Glidden (14)

Peeso Reamers (14)

Paste Carrier (15)

Nickel Titanium Instruments (17-21)

K - Files (18-19)

Super Files III (20-21)

Obturation (23-30)

Gutta Percha (24-27)

Paper Points (28-30)

Accessories (31-38)

Calcium Hydroxide Paste (32-33)

Concord Aspirator (34)

Roeko Paratis Cotton Rolls (35)

Easy Pneumonic Crown & Bridge Remover (36)

Universal Post Remover (38)

Romi Apex Locator (39)

Dentatus (39-68)

Post Systems (39-56)

Surtex Posts (41-47)

Luminex Posts (48-49)

Luscent Anchors (50-51)

Twin Luscent Anchors (52-53)

Reamers (54-55)

Post Keys (55)

Profin (56-68)

Hand Instrument (58)

Contra Angles & Kits (59-61)

Tips (62-68)

Terms and Conditions



Prices: Subject to change without notice.

Shipping: Merchandise shipped UPS.

Claims: Lost or damaged Merchandise must be filed immediately with the delivery company (UPS) as they are responsible for goods delivered to them by Schwed. Claims for shortages, incorrect billing, or delivery of incorrect items must be made to Schwed within 3 days of receipt of the package.

Changes: Terms and Conditions are subject to change at any time.

Return Policy

Return authorization required. No return will be accepted after 30 days from the date of invoice. Original invoice number must accompany all returns.

Merchandise must be returned in its original unopened packaging. Return costs must be prepaid by the customer. There will be a 15% (\$5 minimum) Restocking Charge on all returns. Credit will not be applied to Merchandise not meeting these conditions.

Stainless Steel Instruments



K - Files (4-5)
Hedstrom Files (6-7)
Reamers (8-9)
C - Files (10)
Flexicut Files (11)
Barbed Broaches (12)
Finger Pluggers (13)
Finger Spreaders (13)
Gates Glidden (14)
Peeso Reamers (14)
Paste Filler (15)

K-Files

Schwed K Files have more than two times as many flutes per millimeter as reamers, and the flutes have a greater cutting angle. Because of this design, K files are more efficient instruments for the enlarging of the canal for the removal of dentine and debris. K Files can be used either in a filing motion or in a quarter-turn filing motion.

Like reamers, they should be rotated only in a clockwise direction and should be discarded if any distortion of the flute pattern becomes visible. They must be used in strict sequential order.

Instrument Specifications:



- Blades are made of a special stainless steel alloy, which is endowed with a high degree of elasticity and tensile strength
- All handles are made of plastic, which are ISO color-coded and numbered according to size
- Pre-inserted directional endo stops are color-coded to indicate length (19mm Blue, 21mm White, 25mm Yellow, 28mm Rust, 31mm Black)
- Each instrument is designed for single use
- Made in Germany
- 6 pcs per pkg

	Length	Available Sizes
	19mm	06, 08, 10, 15, 20, 25, 30, 35, 40, 15-40asst, 45
	21mm	06, 08, 10, 12.5, 15, 17.5, 20, 25, 30, 35, 40, 15-40asst, 45, 50, 55, 60, 70, 80, 45-80asst, 90, 100, 110, 90-110asst
	25mm	06, 08, 10, 12.5, 15, 17.5, 20, 25, 30, 35, 40, 15-40asst, 45, 50, 55, 60, 70, 80, 45-80asst, 90, 100, 110, 90-110asst, 120, 130, 140, 120-140asst
	28mm	06, 08, 10, 15, 20, 25, 30, 35, 40, 15-40asst, 45, 50, 55, 60, 70, 80, 45-80asst, 90, 100, 110, 90-110asst
	31mm	06, 08, 10, 15, 20, 25, 30, 35, 40, 15-40asst, 45, 50, 55, 60, 70, 80, 45-80asst, 90, 100, 110, 90-110asst, 120, 130, 140, 120-140asst

K-File Item Numbers

	Size	19mm	21mm	25mm	28mm	31mm
	06	631906	632106	632506	632806	633106
	08	631908	632108	632508	632808	633108
	10	631910	632110	632510	632810	633110
	12.5		632112.5	632512.5		
	15	631915	632115	632515	632815	633115
	17.5		6321175	6325175		
	20	631920	632120	632520	632820	633120
	25	631925	632125	632525	632825	633125
	30	631930	632130	632530	632830	633130
	35	631935	632135	632535	632835	633135
	40	631940	632140	632540	632840	633140
	45	631945	632145	632545	632845	633145
	50		632150	632550	632850	633150
	55		632155	632555	632855	633155
	60		632160	632560	632860	633160
	70		632170	632570	632870	633170
	80		632180	632580	632880	633180
	90		632190	632590	632890	633190
	100		6321100	6325100	6328100	6331100
	110		6321110	6325110	6328110	6331110
	120			6325120		6331120
	130			6325130		6331130
	140			6325140		6331140
	15-40	63191540	63211540	63251540	63281540	63311540
	45-80		63214580	63254580	63284580	63314580
	90-110		632190110	632590110	632890110	633190110
	120-140			6325120140		6331120140

Hedstrom Files

Schwed Hedstrom Files have one continuous cutting edge which runs spirally from the tip of the blade toward the handle. This edge is extremely sharp which makes the Hedstrom File a very effective cutting instrument. The dentine scraped from the canal walls is collected in the chip space on the blade and is efficiently removed. They do not cut with their tips and have no cutting action on insertion. They are effective only on the pull stroke. They must be used in strict sequential order. Care must be taken to avoid rotating the instrument in order to prevent the possibility of a corkscrew effect, or instrument distortion.

Instrument Specifications:

- Blades are made of a special stainless steel alloy which is endowed with a high degree of elasticity and tensile strength
- All handles are made of plastic which are ISO color-coded and numbered according to use
- Pre-inserted directional endo stops are color-coded to indicate length (19mm Blue, 21mm White, 25mm Yellow, 28mm Rust, 31mm Black)
- Each instrument is designed for single use
- Made in Germany
- 6 pcs per pkg



	Length	Available Sizes
	21mm	08, 10, 15, 20, 25, 30, 35, 40, 15-40asst, 45, 50, 55, 60, 70, 80, 45-80asst, 90, 100, 110, 90-110asst, 120, 130, 140
	25mm	08, 10, 15, 20, 25, 30, 35, 40, 15-40asst, 45, 50, 55, 60, 70, 80, 45-80asst, 90, 100, 110, 90-110asst, 120, 130, 140, 120-140asst
	28mm	08, 10, 15, 20, 25, 30, 35, 40, 15-40asst, 45, 50, 55, 60, 70, 80, 45-80asst, 90, 100, 110, 90-110asst
	31mm	08, 10, 15, 20, 25, 30, 35, 40, 15-40asst, 45, 50, 55, 60, 70, 80, 45-80asst, 90, 100, 110, 90-110asst, 120, 130, 140, 120-140asst

Hedstrom File Item Numbers

	Size	21mm	25mm	28mm	31mm
	08	732108	732508	732808	733108
	10	732110	732510	732810	733110
	15	732115	732515	732815	733115
	20	732120	732520	732820	733120
	25	732125	732525	732825	733125
	30	732130	732530	732830	733130
	35	732135	732535	732835	733135
	40	732140	732540	732840	733140
	45	732145	732545	732845	733145
	50	732150	732550	732850	733150
	55	732155	732555	732855	733155
	60	732160	732560	732860	733160
	70	732170	732570	732870	733170
	80	732180	732580	732880	733180
	90	732190	732590	732890	733190
	100	7321100	7325100	7328100	7331100
	110	7321110	7325110	7328110	7331110
	120	7321120	7325120		7331120
	130	7321130	7325130		7331130
	140	7321140	7325140		7331140
	15-40	73211540	73251540	73281540	73311540
	45-80	73214580	73254580	73284580	73314580
	90-110	732190110	73259010	732890110	733190110
	120-140		7325120140		7331120140

Reamers

Schwed Reamers have their cutting edge along the extreme edge of their flutes.

This instrument is used to shave dentine from the canal wall. The spiral flutes collect the dentinal shavings and carry them away from the apex. They are used by rotating the instrument in a clockwise direction during insertion. They must be used in strict sequential order. If the reamer binds, it should be backed-off and reinserted until it can rotate freely and easily. Any distortion of the flute pattern is an indication to discard that instrument. Reamers have no cutting effect as they are withdrawn from the canal.



Instrument Specifications:

- Blades are made of a special stainless steel alloy which is endowed with a high degree of elasticity and tensile strength
- All handles are made of plastic, which are ISO color-coded and numbered according to size
- Pre-inserted directional endo stops are color-coded to indicate length (21mm White, 25mm Yellow, 28mm Rust, 31mm Black)
- Each instrument is designed for single use
- Made in Germany
- 6 pcs per pkg

	Length	Available Sizes
	21mm	06, 08, 10, 15, 20, 25, 30, 35, 40, 15-40asst, 45, 50, 55, 60, 70, 80, 45-80asst, 90, 100, 110, 90-110asst
	25mm	06, 08, 10, 15, 20, 25, 30, 35, 40, 15-40asst, 45, 50, 55, 60, 70, 80, 45-80asst, 100, 90-110asst, 120-140asst
	28mm	15, 20, 25, 30, 35, 40, 15-40asst., 45, 50, 55, 60, 70, 80, 45-80asst, 80, 90, 100, 110, 90-110asst, 130, 140
	31mm	06, 08, 10, 15, 20, 25, 30, 35, 40, 15-40asst, 45, 50, 55, 60, 70, 80, 90, 100, 110, 120, 130, 140

Reamer Item Numbers

	Size	21mm	25mm	28mm	31mm
	06	532106	532506		533106
	08	532108	532508		533108
	10	532110	532510		533110
	15	532115	532515	532815	533115
	20	532120	532520	532820	533120
	25	532125	532525	532825	533125
	30	532130	532530	532830	533130
	35	532135	532535	532835	533135
	40	532140	532540	532840	533140
	45	532145	532545	532845	533145
	50	532150	532550	532850	533150
	55	532155	532555	532855	533155
	60	532160	532560	532860	533160
	70	532170	532570	532870	533170
	80	532180	532580	532880	533180
	90	532190		532890	533190
	100	5321100	5325100	5328100	5331100
	110	5321110		5328110	5331110
	120			5328120	5331120
	130			5328130	5331130
	140			5328140	5331140
	15-40	53211540	53251540	53281540	53311540
	45-80	53214580	53254580	53284580	53314580
	90-110	532190110	532590110	532890110	533190110
	120-140		5325120140		5331120140

C-Files

Schwed “C” files are specially designed for extremely curved or calcified canals. Their stainless steel shafts are heat treated to give them extra body. Their steel blade has extra firmness where the stiff size #06 has the feel of a #10.

Instrument Specifications:

- Blades are made of specially treated high-grade stainless steel
- All handles are made of plastic, which are ISO color-coded and have “C” imprinted on the top of each handle
- Pre-inserted directional endo stops are color-coded to indicate length (21mm White, 25mm Yellow)
- Each instrument is designed for single use
- Made in Germany
- 6 pcs per pkg



	Length	Available Sizes
	21mm	06, 08, 10, 12.5, 15
	25mm	06, 08, 10, 12.5, 15

C-File Item Numbers

	Size	21mm	25mm
	06	652106	652506
	08	652108	652508
	10	652110	652510
	12.5	652112.5	652512.5
	15	652115	652515

Flexicut Files

Schwed Flexicut Files have a triangular cross section instead of the normal square one. This creates a smaller cross section for maximum flexibility. The triangular shape results in 3 sharp cutting edges with an acute rake angle for optimal cutting efficiency and a larger chip space to remove dentinal debris more efficiently. The non aggressive tip avoids ledging, even in the narrowest most tortuous canals. Flexicut will follow any canal curvature all the way to the apex and will clean the entire pulp canal.

Instrument Specifications:

- Blades are made of a special stainless steel alloy, which is endowed with a high degree of elasticity and tensile strength
- All handles are made of plastic, which are ISO color-coded with “F” imprinted on the top
- Pre-inserted directional endo stops are color-coded to indicate length (21mm White, 25mm Yellow)
- Each instrument is designed for single use
- Made in Germany
- 6 pcs per pkg



	Length	Available Sizes
	21mm	15, 20, 25, 30, 35, 40, 15-40asst
	25mm	15, 20, 25, 30, 35, 40, 15-40asst

Flexicut File Item Numbers

	Size	21mm	25mm
	15	642115	642515
	20	642120	642520
	25	642125	642525
	30	642130	642530
	35	642135	642535
	40	642140	642540
	15-40	64211540	64251540

Barbed Broaches

Schwed Barbed Broaches have a thin design that makes them flexible and efficient at removing pulp from the canal. These Broaches have barbs along the shaft. They have plastic handles that are ISO color coded to indicate size. 10 per card.



Instrument Specifications:

- Blades are made of high-grade stainless steel
- ISO color-coded handles
- Length of 21mm
- Each instrument is designed for single use
- Made in China
- 10 pcs per pkg

	Size	Item No.
	XXXX-F (#0)	333PXXXXF
	XXX-F (#1)	333PXXXF
	XX-F (#2)	333PXXF
	X-F (#3)	333PXF
	F (#4)	333PF
	M (#5)	333PM
	C (#6)	333PC
	Assorted	333PA

Assorted Pack contains:

1 - C
 1 - M
 2 - F
 2 - XF
 2 - XXF
 1 - XXXF
 1 - XXXXF

Finger Pluggers

Schwed Finger Pluggers have a smooth conical shape with a flat tip which is designed for both vertical and lateral compaction of Gutta Percha or other filling materials. Plastic handles that are ISO color coded to indicate size. Made in China 6 per pkg.

Finger Plugger Item Numbers

	Size	21mm	25mm	31mm
	15	992115	992515	993115
	20	992120	992520	993120
	25	992125	992525	993125
	30	992130	992530	993130
	35	992135	992535	993135
	40	992140	992540	993140
	15-40	99211540	99251540	99311540



** Special ordersizes 06, 08, 10, 45, 50, 55, 60, 70, 80 **

Finger Spreaders

Schwed Finger Spreaders have a smooth conical shape with a pointed tip for compaction of Gutta Percha or other filling materials. Plastic handles that are ISO color coded to indicate size. Made in China. 6 per pkg.

Finger Spreader Item Numbers

	Size	21mm	25mm	31mm
	15	952115	952515	953115
	20	952120	952520	953120
	25	952125	952525	953125
	30	952130	952530	953130
	35	952135	952535	953135
	40	952140	952540	953140
	15-40	95211540	95251540	95311540



** Special ordersizes 06, 08, 10, 45, 50, 55, 60, 70, 80 **

Gates Glidden Drills

Schwed Gates Glidden Drills are side cutting engine instruments that can be used either as an orifice cleaner or to shape the coronal half of the canal. The number of bands on the handle represent the size of the instrument. Made from Stainless Steel to allow for flexibility in order to prevent breakage. Made in Germany 6 per pkg.

Gates Glidden Item Numbers

Size	28mm	32mm
1	91S1	911
2	91S2	912
3	91S3	913
4	91S4	914
5	91S5	915
6	91S6	916
1-6	91SA	91A



Peeso Reamers

Schwed Peeso Reamers are side cutting engine instruments that are used to funnel the coronal half of the canal. They are also helpful in preparing space for a dowel or post. Size is indicated by the number of bands on the handle of the instrument. They are made from Stainless Steel to allow for flexibility to prevent breakage. Made in Germany. 6 per pkg.

Peeso Reamer Item Numbers

Size	28mm	32mm
1	47S1	471
2	47S2	472
3	47S3	473
4	47S4	474
5	47S5	475
6	47S6	476
1-6	47SA	47A



Paste Carriers

Schwed Paste Carriers are made up of a fine wire that is twisted into the shape of a reverse spiral. They carry cement or filling paste apically as they rotate in a clockwise direction. They can be used either with a slow-speed contra-angle, or as a hand instrument.

Instrument Specifications:

- Made to ISO standards
- Blade is made of a high grade stainless steel
- Handles are made of plastic and are color coded to indicate size
- Made in China
- 4 pieces per package



Paste Carrier Item Numbers

	Size	21mm	25mm	29mm
	25	932125	932525	932925
	30	932130	932530	932930
	35	932135	932535	932935
	40	932140	932540	932940
	15-40	93211540	93251540	93291540

Nickel Titanium Instruments



K - Files (18-19)
Super Files III (20-21)

NiTi K-Files

With the Schwed Niti K File you can enlarge just about any canal curvature more safely, more accurately and faster with more flexibility and stress resistance.



Instrument Specifications:

- Blades are made of Nickel Titanium (NiTi) with exceptional flexibility.
- All handles are plastic, which are ISO color coded and numbered according to size.
- Pre-inserted directional endo stops are color-coded to indicate length (21mm White, 25mm Yellow)
- .02 Taper
- Made in China
- 6 pcs per pkg

Compared with Conventional Files:

- Greater the flexibility to follow even extreme canal curvatures.
- Increased stress resistance to prevent breakage.
- Lasts longer
- 100% memory. Pops back into linear alignment immediately without distortion.
- No tendency to straighten the canal. Follows the existing canal path.
- Does not change shape of the canal during instrumentation.
- Does not transpose the apex or cause ledging.

	Length	Available Sizes
	21mm	08, 10, 15, 20, 25, 30, 35, 40, 15-40asst, 45, 50, 55, 60, 70, 80, 45-80asst
	25mm	08, 10, 15, 20, 25, 30, 35, 40, 15-40asst, 45, 50, 55, 60, 70, 80, 45-80asst
	28mm	10, 15, 20, 25, 30, 35, 40, 15-40asst, 45, 50, 55, 60, 70, 80, 45-80asst
	31mm	10, 15, 20, 25, 30, 35, 40, 15-40asst, 45, 50, 55, 60, 70, 80, 45-80asst

NiTi K-File Item Numbers

	Size	21mm	25mm	28mm	31mm
	08	NK2108	NK2508		
	10	NK2110	NK2510	NK2810	NK3110
	15	NK2115	NK2515	NK2815	NK3115
	20	NK2120	NK2520	NK2820	NK3120
	25	NK2125	NK2525	NK2825	NK3125
	30	NK2130	NK2530	NK2830	NK3130
	35	NK2135	NK2535	NK2835	NK2135
	40	NK2140	NK2540	NK2840	NK3140
	45	NK2145	NK2545	NK2845	NK3145
	50	NK2150	NK2550	NK2850	NK3150
	55	NK2155	NK2555	NK2855	NK3155
	60	NK2160	NK2560	NK2860	NK3160
	70	NK2170	NK2570	NK2870	NK3170
	80	NK2180	NK2580	NK2880	NK3180
	15-40	NK211540	NK251540	NK281540	NK311540
	45-80	NK214580	NK254580	NK284580	NK314580



Super Files III

Super Files III are made of high quality Nickel Titanium. Their design makes them highly flexible for moderate to severe bending of the root canal for use on curved, narrow, calcified, and other complex root canals. The secure round tip (0.4mm) avoids lateral perforation or ledging. Minimal risk design of instrument separation and step side penetration. Match these files with the Protaper Gutta Percha.



Instrument Specifications:

- Made of Nickel Titanium
- High Quality Controlled Memory Wire
- Heat Treated (Blue Color)
- Made in China
- 6 pcs per box or blister



Operation Process		
PROBE	#10	
	#15	
OPENING	SX	
SHAPING	S1	
	S2	
FINISHING	F1	
	F2	
	F3	
ENLARGING	F4	
	F5	

Super Files III Item Numbers

	Size	Taper (%)	Tip (mm)	19mm	21mm	25mm	31mm
	SX	3.5-19	0.19	NRS319SX			
	S1	2-11	0.17		NRS321S1	NRS325S1	NRS331S1
	S2	4-11.5	0.20		NRS321S2	NRS325S2	NRS331S2
	F1	7	0.20		NRS321F1	NRS325F1	NRS331F1
	F2	8	0.25		NRS321F2	NRS325F2	NRS331F2
	F3	9	0.30		NRS321F3	NRS325F3	NRS331F3
	SX-F3 Asst				NRS321AS	NRS325AS	NRS331AS
	F4	6	0.40		NRS321F4	NRS325F4	NRS331F4
	F5	5	0.50		NRS321F5	NRS325F5	NRS331F5

Matching Protaper Gutta Percha can be found on page 26

Matching Protaper Paper Points can be found on page 30

Obturation



Gutta Percha Points (24-25)
Accessory Gutta Percha (26)
Protaper Gutta Percha (26)
Gutta Percha Pellets (27)
Paper Points (28-29)
Accessory Paper Points (30)
Protaper Paper Points (30)

Gutta Percha Points

Schwed Meta Gutta Percha Points come in a sliding package for easy dispensing and selection. They are Cadmium free and color coded with non-toxic pigments.



Product Specifications:

- ISO .02, .04, and .06 taper
- Sizes 10 through 140
- **60** or **120** pieces per package
- Stiff and Flexible
- Highly Radiopaque
- Precision hand rolled
- Suitable for various techniques
- Made in Korea

Taper	Available Sizes
.02	10, 15, 20, 25, 30, 35, 40, 15-40asst, 45, 50, 55, 60, 70, 80, 45-80asst, 90, 100, 110, 120, 130, 140
.04	15, 20, 25, 30, 35, 40, 15-40asst, 45, 50, 55, 60, 70, 80, 45-80asst
.06	15, 20, 25, 30, 35, 40, 15-40asst, 45, 50, 55, 60, 70, 80, 45-80asst

Gutta Percha Point Item Numbers

	Size	.02 Taper	.04 Taper	.06 Taper
	10	GP2-10		
	15	GP2-15	GP4-15	GP6-15
	20	GP2-20	GP4-20	GP6-20
	25	GP2-25	GP4-25	GP6-25
	30	GP2-30	GP4-30	GP6-30
	35	GP2-35	GP4-35	GP6-35
	40	GP2-40	GP4-40	GP6-40
	45	GP2-45	GP4-45	GP6-45
	50	GP2-50	GP4-50	GP6-50
	55	GP2-55	GP4-55	GP6-55
	60	GP2-60	GP4-60	GP6-60
	70	GP2-70	GP4-70	GP6-70
	80	GP2-80	GP4-80	GP6-80
	90	GP2-90		
	100	GP2-100		
	110	GP2-110		
	120	GP2-120		
	130	GP2-130		
	140	GP2-140		
	15-40	GP2-1540	GP4-1540	GP6-1540
	45-80	GP2-4580	GP4-4580	GP6-4580

Accessory Gutta Percha



Product Specifications:

- Cadmium Free
- Sizes XF through ML
- **60** or **120** pieces per package
- Stiff and Flexible
- Highly Radiopaque
- Precision hand rolled
- Made in Korea

Size	XF	FF	MF	F	FM	M	ML
Item No.	GP-XF	GP-FF	GP-MF	GP-F	GP-FM	GP-M	GP-ML

Protaper Gutta Percha



Product Specifications:

- Sizes F1 - F3
- **60** pieces per package
- Stiff and Flexible
- Highly Radiopaque
- Precision hand rolled
- Length markings
- Made in Korea

Size	F1	F2	F3	F1-F3
Item No.	GP-F1	GP-F2	GP-F3	GP-F1F3

Gutta Percha Pellets

Schwed Meta Gutta Percha Pellets are uniformly rolled from the finest gutta percha to extremely close tolerances. They are used for obturation

Product Specifications:

- Natural
- Maximum Flowability
- Premium Melting Temp. at 160-200 degrees Celsius
- Needle Size: 23 or 25 gauge
- 100 Pellets per pkg
- Made in Korea



Item Number: GP-OB



Paper Points

Schwed Meta Absorbent Paper Points are primarily used to dry the root canal. They can also be helpful in introducing various medications and sealers into the canal.

Sterile points can be used to culture the canal.



Product Specifications:

- ISO .02, .04, and .06 taper in sizes 10 through 140
- 60, 100, or 200 pieces per package
- Made of highly absorbent paper.
- Tightly hand rolled ensuring better retention of size and shape when wet
- Made in Korea

Taper	Available Sizes
.02	10, 15, 20, 25, 30, 35, 40, 15-40asst, 45, 50, 55, 60, 70, 80, 45-80asst, 90, 100, 110, 120, 130, 140
.04	15, 20, 25, 30, 35, 40, 15-40asst, 45, 50, 55, 60, 70, 80, 45-80asst
.06	15, 20, 25, 30, 35, 40, 15-40asst, 45, 50, 55, 60, 70, 80, 45-80asst

Paper Point Item Numbers

	Size	.02 Taper	.04 Taper	.06 Taper
	10	PP2-10		
	15	PP2-15	PP4-15	PP6-15
	20	PP2-20	PP4-20	PP6-20
	25	PP2-25	PP4-25	PP6-25
	30	PP2-30	PP4-30	PP6-30
	35	PP2-35	PP4-35	PP6-35
	40	PP2-40	PP4-40	PP6-40
	45	PP2-45	PP4-45	PP6-45
	50	PP2-50	PP4-50	PP6-50
	55	PP2-55	PP4-55	PP6-55
	60	PP2-60	PP4-60	PP6-60
	70	PP2-70	PP4-70	PP6-70
	80	PP2-80	PP4-80	PP6-80
	90	PP2-90		
	100	PP2-100		
	110	PP2-110		
	120	PP2-120		
	130	PP2-130		
	140	PP2-140		
	15-40	PP2-1540	PP4-1540	PP6-1540
	45-80	PP2-4580	PP4-4580	PP6-4580

Accessory Paper Points



Product Specifications:

- Accessory sizes XF through XC
- **100** or **200** pieces per package
- Made of highly absorbent paper
- Tightly hand rolled ensuring better retention of size and shape when wet
- Made in Korea

Size	XXF	XF	F	M	C	XC
Item No.	PP-XXF	PP-XF	PP-FF	PP-M	PP-C	PP-XC

Protaper Paper Points



Product Specifications:

- Protaper F1-F3
- **60** pieces per package
- Made of highly absorbent paper
- Length markings
- Tightly hand rolled ensuring better retention of size and shape when wet
- Made in Korea

Size	F1	F2	F3	F1-F3
Item No.	PP-F1	PP-F2	PP-F3	PP-F1F3

Accessories



Calcium Hydroxide Paste (32-33)
Concord Aspirator (34)
Roeko Paratis Cotton Rolls (35)
Easy Pneumatic Crown & Bridge Remover (36)
Universal Post Remover (37)
RomiApex Locator A-15 (38)

Metapaste Plus

Metapaste Plus is a Calcium Hydroxide Paste that comes in premixed tubes for quick and easy treatment with minimal clean up. Used for applications of weeping canals, apexification, temporary root canal filling, and internal and external root resorption. Made in Korea.



Item Number: CHP-01

Attributes:

- *Components:* 2g syringe 1ea, Disposable Tips 20ea
- *Color:* White
- *Composition:* Calcium Hydroxide, Barium Sulfate, Purified Water
- *Shelf life:* Two years from production date

Advantages:

- Faster medicinal effect and easily removable
- PPG Free to reduce post-op sensitivity
- Antibacterial with high pH: 12.15 (1 Day) / 12.35 (14 Days)
- Improved Flowability enhances convenience and accuracy
- High Radiopacity of more than 3Almm

Metapex Plus

Metapex Plus is a Calcium Hydroxide Paste that comes in premixed tubes for quick and easy treatment with minimal clean up. Used for applications of weeping canals, apexification, temporary root canal filling, and internal and external root resorption. Made in Korea.



Item Number: CHP-02

Attributes:

- *Components:* 2.2g syringe 1ea, Disposable Tips 20ea
- *Color:* Yellow
- *Composition:* Calcium Hydroxide, Iodoform, Silicon Oil
- *Shelf life:* Two years from production date

Advantages:

- Antibacterial with high pH: 12.00 (1 Day) / 12.40 (14 Days)
- Improved Flowability enhances convenience and accuracy
- High Radiopacity of 5.4Almm

Concord Aspirator

The Schwed Concorde Aspirator has a unique, innovative design advance offering numerous benefits compared to conventional aspirators. The 'High-Tech' flange design and bias angled orifice prevents 'suction-grab' while the larger opening removed debris and liquids faster and quieter. The overall size and shape also make Concorde a perfect retractor to become a dual functioning instrument. Slightly angled, the 5" length fits comfortably into the hand for versatile manipulation. It is autoclavable to 130 degrees C, fits all standard aspirator systems and has a non-glare matte finish.



Product Specifications:

- Long and fits most standard aspirator systems
- Can be sterilized in autoclave or chemiclave up to 130 degrees C.
- Aerodynamically Designed to Prevent 'Suction-Grab' of Soft Tissue
- Dual Functioning Instrument
- Used as a combined aspirator and retractor
- Low noise suction
- Made in Germany
- 10 aspirators to a pack



Size	Item No.
11mm	450211
16mm	450216

Roeko Cotton Rolls

Schwed Parotis Cotton Rolls provide an extensive, dry operative field. The Parotis Roll’s spring-action occludes the parotid duct, thereby stopping the saliva flow. This virtually eliminates the need for saliva ejectors, simultaneously keeping the patient’s mouth open and retracting the patient’s cheek.

Product Specifications:

- Each roll consists of a spring completely encased in a highly absorbent cotton
- Pharmaceutical quality raw materials. Tested for purity, foreign fibers, absorbency, and surfactant residues
- Easy to insert and comfortable for the patient
- Disposable
- Cotton designed to not unravel during use
- Available in four sizes
- Made in Germany
- 100 pieces per package



Size	Type	Length	Diameter	Item No.
1	Child	80mm	9mm	160001
2	Child	80mm	10mm	160002
3	Adult	100mm	10mm	160003
5	Adult	100mm	8mm	160005

Easy Pneumatic Crown & Bridge Remover



Item Number: 0100

Kit Includes:


- Crown & Bridge Remover
- Bridge Remover Tip
- Single Crown Remover Tip
- Single Crown Remover Contrangled Tip
- Open end Adjustment wrench
- Flexi-wire Kit
- Instructional video / CD ROM
- Storage / Carrying case

Advantages:

- Less chance of tooth fracture than traditional crown and bridge removers
- Restorations removed intact and reusable
- No porcelain damage
- Significantly reduced chair time
- Easy fingertip control
- Maintenance free
- Uses a controlled, predictable force
- Uses a compressed air with adjustable speed force
- High patient acceptance
- Requires no anesthesia
- No need for an assistant
- Connect directly to 4/5 or 2/3-hole handpiece connector
- OSHA compliant: fully autoclavable

Universal Post Remover

The new improved Gonon Universal Post Remover enables the practitioner to remove whole or broken posts from root canals prior to retreatment. It can be used on posts made of steel, gold and stainless steel. The kit is made of stainless steel and comes with four different tool sizes.



Kit Includes:

- 1 Pointer Drill
- 8 Trephine Burs
- 8 CCW Tubular Taps
- 1 Plier with Screwing Knob
- 2 Sets of Washers
- 2 Supporting Tubes
- 1 Aluminum Box

Advantages:

- For whole or broken posts
- Use on posts made of steel, gold, gold-plated or stainless steel

Features:

- Made of Stainless Steel
- 4 Different Sizes
- 3 Easy steps to post removal

	Description	Item No.
	Universal Post Remover Kit	T-POST-UPR
	Trephine Bur for Post Remover #1	T-BUR-1
	Trephine Bur for Post Remover #2	T-BUR-2
	Trephine Bur for Post Remover #3	T-BUR-3
	CCW Tap #1 for Screw Post	T-CCW-TAP-1
	CCW Tap #2 for Screw Post	T-CCW-TAP-2
	CCW Tap #3 for Screw Post	T-CCW-TAP-3
	CCW Tap #4 for Screw Post	T-CCW-TAP-4

RomiApex Locator A-15

Apex Locator has advanced technology and a user friendly design.

Basic Kit Includes:

- RomiApex A-15 locator unit
- Measuring Cable (1)
- File Clips (2)
- Lip Clips (5)
- Touch Probe (1)



Product Advantages:

- Advanced Technology
- Compact size
- User Friendly design
- Fully Automatic Measurements
- High precision of apical foramen localization
- Works in dry or wet canals
- Audio feedback with volume control
- Automatic turn-off

Displayed Numerical Values-Distance to Apex

2.0mm, 1.5mm, 1.0mm, 0.75mm, 0.50mm, 0.25mm, 0.00mm, FO (over Foramen)

Item Number: Apex A-15

Specifications:

- Dimensions: 55mm Width X 92mm High X 16mm Thickness
- Weight: 100gr
- Display Type: Color LCD
- Display Size: 51mm X 38mm
- Power Supply: 1 AAA Alkaline Batteries

Dentatus Post Systems




Surtex Posts (41-47)
Luminex Posts (48-49)
Luscent Anchors (50-51)
Twin Luscent Anchors (52-53)
Reamers (54-55)
Post Keys (55)

Dentatus Surtex Posts

Surtex® Classic Posts are reliable solutions for post and core build-ups. The high-quality, prefabricated posts are dependable and economical.

Available in gold-plated, titanium, and stainless steel materials, they are intended for passive cementation in root canals. The thread increases retention with these posts as well as increased surface area and easy conal re-access if needed, and the innovative surface texture delivers 200% more mechanical retention as compared to non-textured posts. This textured surface not only minimizes the risk of posts loosening, but also masks post exposure through restorative materials. The crosscut head can be expanded with the cross cut key to provide mechanical retention to secure the core build-up. The anatomical canal adaptation parallel sided with a tapered apical end, minimizes the risk of root fracture.

The Schwed Dentatus fully integrated One-Stop Classic Post System is your all in one kit for posts. The One-Stop Post Systems come with a variety of posts, accessories, and reamers to match.



Gold Plated

Surtex Gold Plated (RST) posts are recommended when an economical alternative is required.

Gold Plated Post Kit, 18 Posts

Item No. RST-18



Includes:

- 18 Gold-Plated Posts
 - 1 of each
 - Short #1, #2
 - Medium #1, #4, #5, #6
 - Long #1, #2, #3, #4, #5, #6
 - Extra Long #4, #6
 - 2 of each
 - Medium #2, #3
- 6 Reamers
 - 1 of each
 - #1, #2, #3, #4, #5, #6
- 1 Hollow Key, 1 Cross Key, 20 Length Stops

Gold Plated Post Kit, 60 Posts

Item No. RST-60



Includes:

- 60 Gold-Plated Posts
 - 3 of each
 - Short #1, #2, #3, #4, #5, #6
 - Medium #1, #2, #3, #4, #5, #6
 - Long #1, #2, #3, #4, #5, #6
 - Extra Long #4, #6
- 6 Reamers
 - 1 of each
 - #1, #2, #3, #4, #5, #6
- 1 Hollow Key, 1 Cross Key, 20 Length Stops

Gold Plated

Gold Plated Post Refills

Refills contain 12 posts per pkg



Gold Plated Post Refill Item Numbers

Size	Diameter	Short 7.8mm	Medium 9.3mm	Long 11.8mm	Extra Long 14.2mm	Super Long 17.0mm
1	1.05mm	RST-S1	RST-M1	RST-L1		
2	1.20mm	RST-S2	RST-M2	RST-L2		
3	1.35mm	RST-S3	RST-M3	RST-L3	RST-EL3	RST-SL3
4	1.50mm	RST-S4	RST-M4	RST-L4	RST-EL4	RST-SL4
5	1.65mm	RST-S5	RST-M5	RST-L5	RST-EL5	RST-SL5
6	1.80mm	RST-S6	RST-M6	RST-L6	RST-EL6	RST-SL6

Titanium

Surtex Titanium (TST) posts are recommended for optimal biocompatibility with surrounding tissue.

Titanium Post Kit, 60 Posts Item No. TST-60



Includes:

- 60 Titanium Posts
 - 3 of each
 - Short #1, #2, #3, #4, #5, #6
 - Medium #1, #2, #3, #4, #5, #6
 - Long #1, #2, #3, #4, #5, #6
 - Extra Long #4, #6
- 6 Reamers
 - 1 of each
 - #1, #2, #3, #4, #5, #6
- 1 Hollow Key, 1 Cross Key, 20 Length Stops

Titanium

Titanium Post Refills

Refills contain 12 posts per pkg




Titanium Post Refill Item Numbers

Size	Diameter	Short 7.8mm	Medium 9.3mm	Long 11.8mm	Extra Long 14.2mm	Super Long 17.0mm
1	1.05mm	TST-S1	TST-M1	TST-L1		
2	1.20mm	TST-S2	TST-M2	TST-L2		
3	1.35mm	TST-S3	TST-M3	TST-L3	TST-EL3	TST-SL3
4	1.50mm	TST-S4	TST-M4	TST-L4	TST-EL4	TST-SL4
5	1.65mm	TST-S5	TST-M5	TST-L5	TST-EL5	TST-SL5
6	1.80mm	TST-S6	TST-M6	TST-L6	TST-EL6	TST-SL6

Stainless Steel

Surtex Stainless Steel (SST) posts are recommended for all retentive purposes where extra strength is required.

Stainless Steel Post Kit, 60 Posts Item No. SST-60



Includes:

- 60 Stainless Steel Posts
 - 3 of each
 - Short #1, #2, #3, #4
 - Medium #1, #2, #3, #4
 - Long #3, #4, #5, #6
 - Extra Long #3, #4, #5 #6
 - Super Long #3, #4, #5, #6
- 6 Reamers
 - 1 of each
 - #1, #2, #3, #4, #5, #6
- 1 Hollow Key, 1 Cross Key, 20 Length Stops

Stainless Steel

Stainless Steel Post Refills

Refills contain 12 posts per pkg



Stainless Steel Post Refill Item Numbers

Size	Diameter	Short 7.8mm	Medium 9.3mm	Long 11.8mm	Extra Long 14.2mm	Super Long 17.0mm
1	1.05mm	SST-S1	SST-M1	SST-L1		
2	1.20mm	SST-S2	SST-M2	SST-L2		
3	1.35mm	SST-S3	SST-M3	SST-L3	SST-EL3	SST-SL3
4	1.50mm	SST-S4	SST-M4	SST-L4	SST-EL4	SST-SL4
5	1.65mm	SST-S5	SST-M5	SST-L5	SST-EL5	SST-SL5
6	1.80mm	SST-S6	SST-M6	SST-L6	SST-EL6	SST-SL6

Luminex Posts

The Luminex post technique centers and forms a parallel-sided canal for corresponding Dentatus Surtex metal posts. Luminex posts have superior aesthetics for light-cured composite inside the canal that masks metal posts with a reflective tooth colored foundation for modern restorations. Having technique versatility, Luminex smooth and grooved posts may also be used as an impression and castable post pattern in the direct and indirect fabrication of posts. Luminex Posts are available in the Luminex Post System and refills to match.

All too often, fragile, thin-walled teeth present major restorative problems; cast posts or extractions are the only alternative. But now Luminex®, a clear, light-transmitting, plastic post, polymerizes light-cured composites within the entire root canal. After curing, the Luminex post is removed, leaving a ready canal for a corresponding Surtex Post.



Problem

The large defect at coronal root end of the maxillary left central incisor is caused by carious extension into a previously traumatized, immature incisor.



A matching post canal is created in conjunction with composite resin reinforcement of the defective root.




Solution

Light is transmitted into the root canal using light-transmitting plastic post to polymerize the intracanal composite resin.



A prefabricated matching grooved metal post is cemented into the reconstituted post canal.



Luminex Posts


*To be used with *Titanium Surtex Posts* - Refer to page 47

Luminex Post Kit

Item No. LST-T


Includes:

- 30 Luminex Smooth Posts
 - 5 of each
 - #1, #2, #3, #4, #5, #6
- 30 Luminex Grooved Posts
 - 5 of each
 - #1, #2, #3, #4, #5, #6
- 40 Surtex Titanium Posts
 - 5 of each
 - Long #1, #2, #3, #4
 - Super Long #3, #4, #5, #6
- 6 Reamers
 - 1 of each
 - #1, #2, #3, #4, #5, #6
- 1 Hollow Key, 1 Cross Key, 20 Length Stops



Luminex Post Refills

Refills contain 20 posts per pkg




Luminex Post Refill Item Numbers

Size	Diameter	Smooth	Grooved
1	1.05mm	DIP-1	DCP-1
2	1.20mm	DIP-2	DCP-2
3	1.35mm	DIP-3	DCP-3
4	1.50mm	DIP-4	DCP-4
5	1.65mm	DIP-5	DCP-5
6	1.80mm	DIP-6	DCP-6

Luscent Anchors

Luscent® Anchors are made of unidirectional fiberglass rods encapsulated in a strong resin matrix. The rods transmit light, ensuring complete cure of core material. Luscent Anchors also reflect natural hues and eliminate shadowing through composites at the margin, making them ideal for highly esthetic restorations.

In addition to aesthetics, Luscent Anchors feature excellent strength characteristics. Their flexural modulus is approximately 20 GPA, very similar to that of natural dentine. Light or dual cure composites bond to the Luscent Anchors, which allow sufficient time for exact centering and seating of the anchor, as well as a complete fill of the composite inside the canal thus creating a cohesive, strong monobloc foundation. Luscent Anchors can be placed in both narrow and large canals with conventional adhesive techniques and materials.



Luscent Anchors

12 Luscent Anchor Kit

Item No: LUC-K1

Includes:

- 12 Luscent Anchors
 - 4 of each
 - Small, Medium Large
- 3 Reamers
 - 1 of each
 - RUA-4, RUA-5, RUB-6
- 12 Assorted Core Forms-to-Fit



Luscent Anchor Refills

	Item No.	Description
	LUC-XXXS5	5 Luscent Anchors in XXXSmall, 6 Asst. Core Forms-to-Fit
	LUC-XXS5	5 Luscent Anchors in XXSmall, 6 Asst. Core Forms-to-Fit
	LUC-XS5	5 Luscent Anchors in XSmall, 6 Asst. Core Forms-to-Fit
	LUC-S5	5 Luscent Anchors in Small, 6 Asst. Core Forms-to-Fit
	LUC-M5	5 Luscent Anchors in Medium, 6 Asst. Core Forms-to-Fit
	LUC-L5	5 Luscent Anchors in Large, 6 Asst. Core Forms-to-Fit
	LUC-S15	15 Luscent Anchors in Small, 15 Asst. Core Forms-to-Fit
	LUC-M15	15 Luscent Anchors in Medium, 15 Asst. Core Forms-to-Fit
	LUC-L15	15 Luscent Anchors in Large, 15 Asst. Core Forms-to-Fit

Twin Luscent Anchors

The innovative design of the Twin Luscent® Anchors is a visible assurance against accidental debonding of adhesive and resin-core materials. The slim mid-section creates a “physical choke.” The vent groove eliminates air resin entrapment and prevents rotational dislocation. A great combination of light transmission, attractive aesthetics and twice the retention.

They eliminate shadows at the gingival, root and crown interface as well as through thin laminate composite restorations. The Twin Luscent Anchors reflect the surrounding colors and hues, compatible with natural esthetics. The anchor’s cone-shaped end is used in narrower canals, preventing excess removal of dentin or canal wall. The parallel-end is used in wider canals. The parallel canals can be refined with drills, using parallel-canal post techniques. The Anchor’s elasticity (20.1GPA) is within the range of healthy teeth, providing a safe and cohesive resistance to impact.



Twin Luscent Anchors

15 Twin Luscent Anchor Kit

Item No: LUT-K1

Includes:

- 15 Twin Luscent Anchors
 - 5 of each
 - Small, Medium Large
- 3 Reamers
 - 1 of each
 - RUA-4, RUA-5, RUB-6
- 1 Pathfinder
 - RUF-4
- 1 Probos II
 - RUD-4
- 15 Assorted Core Forms-to-Fit



Twin Luscent Anchor Refills

	Item No.	Description
	LUT-XS5	5 Twin Luscent Anchors in XSmall, 6 Asst. Core Forms-to-Fit
	LUT-S5	5 Twin Luscent Anchors in Small, 6 Asst. Core Forms-to-Fit
	LUT-M5	5 Twin Luscent Anchors in Medium, 6 Asst. Core Forms-to-Fit
	LUT-L5	5 Twin Luscent Anchors in Large, 6 Asst. Core Forms-to-Fit
	LUT-S15	15 Twin Luscent Anchors in Small, 15 Asst. Core Forms-to-Fit
	LUT-M15	15 Twin Luscent Anchors in Medium, 15 Asst. Core Forms-to-Fit
	LUT-L15	15 Twin Luscent Anchors in Large, 15 Asst. Core Forms-to-Fit

Standard Classic Reamers

Dentatus Reamers are made in Sweden from hardened stainless steel. They create a perfect match between root canals and posts. Sizes are indicated by the number of grooves on the shank, which correspond to the post size on the Dentatus metal posts.

Standard Classic Reamers are for Surtex, Twin Luscent and Luminex posts.



Standard Classic Reamer Refills

Item No.	Length	Description
RUA	28mm	6 Short Reamers, 1 of each #1 through #6
RUA-1	28mm	3 Short Reamers #1
RUA-2	28mm	3 Short Reamers #2
RUA-3	28mm	3 Short Reamers #3
RUA-4	28mm	3 Short Reamers #4
RUA-5	28mm	3 Short Reamers #5
RUA-6	28mm	3 Short Reamers #6
RUB	33mm	6 Long Reamers, 1 of each #1 through #6
RUB-1	33mm	3 Long Reamers #1
RUB-2	33mm	3 Long Reamers #2
RUB-3	33mm	3 Long Reamers #3
RUB-4	33mm	3 Long Reamers #4
RUB-5	33mm	3 Long Reamers #5
RUB-6	33mm	3 Long Reamers #6

Special Reamers

Probos II Routers

Dentatus Probos II Router improves retention by creating indents in the dentin wall when moved laterally.

Item No.	Description
RUD	4 Probos II Routers, 1 of each #3 through #6
RUD-3	3 Probos II Routers #3
RUD-4	3 Probos II Routers #4
RUD-5	3 Probos II Routers #5
RUD-6	3 Probos II Routers #6



Pathfinder Reamer

Dentatus Pathfinder Reamer is used for removal of root filling material.



Item No.	Description
RUF	6 Pathfinder Reamers, 1 of each #1 through #6
RUF-1	3 Pathfinder Reamers #1
RUF-2	3 Pathfinder Reamers #2
RUF-3	3 Pathfinder Reamers #3
RUF-4	3 Pathfinder Reamers #4
RUF-5	3 Pathfinder Reamers #5
RUF-6	3 Pathfinder Reamers #5

Hollow Key

Item No: NLH



Cross Key

Item No: NLK



Dentatus Profin



Hand Instrument (58)
Directional System (59-68)
Contra Angles & Kits (59-61)
Tips (62-68)

Profin PDH Hand Instrument

The Profin PDH™ hand instrument enables precise finishing work in the most inaccessible areas. It uses Lamineer Tips, which can be set at either end of the instrument. At one end, the tip is perpendicular to the instrument and on the other end the tip is parallel with the instrument. The tip may be set in a locked position, typically for controlled material reduction, or it can be set to freely rotate around its axis to follow the tooth contours.

1 Profin Hand Instrument
w/ 2 Tips

Item No: PDH-SK

Includes:

- 1 Profin PDH Hand Instrument
- 1 LTA-150 Lamineer Diamond Tip
- 1 LAT-75 Lamineer Diamond Tip



Profin Directional System

The Profin® Directional System is the instrument for shaping, contouring and polishing natural dentition and cosmetic restorations with easy access and flexibility. It removes bulk interproximal and subgingival residual bonding materials in all parts of the dental arch. The system comprises a choice of contra-angles and the flat Lamineer® diamond tips, which move in a reciprocating axial direction. With the Lamineer Tips, the range of the system can reach, access, and refinish the entirety of the dental arch.

Ideal for finishing restorations made from diverse materials; reshaping restorations of diverse materials; refining interdental gingival and incisal embrasures; and root planing, scaling, and reshaping of existing restoration margins, extensions, and overhangs.



E-type Profin Contra-Angle for Electric or Air Motors

Item No: PDX



Contra Angles & Kits

E-type Profin Contra-Angle w/ 5 Lamineer Tips

Item No: PDX/SK

Includes:

- 1 E-Type Profin Contra-Angle
- 5 Lamineer Tips
 - 4 Diamond Tips: LTA-150, LTA-S100, LTA-S50, LTA-30
 - 1 Wolfram Tungsten Finishing Tip: LTA-39
- 1 Lamineer Tip Stand w/ Color-Coded Chart
- 20 Plastic Polishing Tips: 10 of EVA-2000, 10 of EVA-5000



E-type Profin Contra-Angle w/ 10 Lamineer Tips

Item No: PDX-K10

Includes:

- 1 E-Type Profin Contra-Angle
- 10 Lamineer Tips
 - 8 Diamond Tips: LTA-S150, LTA-100, LTA-50, LTA-T30, LTA-30, LTA-D30, LTA-S15
 - 2 Wolfram Tungsten Finishing Tips: LTA-S36, LTA-39
- 1 Lamineer Tip Stand w/ Color-Coded Chart
- 100 Plastic Polishing Tips



E-type Profin Contra-Angle w/ 5 IPR Tips

Item No: PDX-P/SK2

Includes:

- 1 E-Type Profin Contra-Angle
- 5 IPR Tips
 - LTA-P5, LTA-P4, LTA-P3, LTA-P2, LTA-P1
- 1 Lamineer Tip Stand w/ Color-Coded Chart
- 20 Plastic Polishing Tips: 10 of EVA-2000, 10 of EVA-5000



Contra Angles & Kits

Doriot Profin Contra-Angle (fits on straight nose cone)

Item No: PDD



Doriot Profin Contra-Angle w/ 10 Lamineer Tips

Item No: PDD-K10

Includes:



- 1 Doriot Profin Contra-Angle
- 10 Lamineer Tips
 - 8 Diamond Tips: LTA-S150, LTA-100, LTA-50, LTA-T30, LTA-30, LTA-D30, LTA-S15
 - 2 Wolfram Tungsten Finishing Tips: LTA-S36, LTA-39
- 1 Lamineer Tip Stand w/ Color-Coded Chart
- 100 Plastic Polishing Tips

Doriot Profin Contra-Angle w/ 5 Lamineer Tips

Item No: PDD-P5


Includes:



- 1 Doriot Profin Contra-Angle
- 4 Lamineer Tips
 - 4 Diamond Tips: LTA-150, LTA-S100, LTA-S50, LTA-30
 - 1 Wolfram Tungsten Finishing Tips: LTA-39
- 1 Lamineer Tip Stand w/ Color-Coded Chart
- 20 Plastic Polishing Tips: 10 of EVA-2000, 10 of EVA-5000

Profin Tips

Lamineer Tips



Lamineer® Tips are offered in different shapes and grits for many applications unique to the Profin® Directional System. They are safe-sided and do not harm soft (gingival) tissues and are available in a wide range of color-coded diamond grit sizes, from coarse (150 microns) for gross removal of all restorative materials, to extra-fine (15 microns) for polishing. These flat, thin tips allow optimal access with the diamond coating only on one side for safety and control.

Wolfram tungsten tips are also available for even finer polishing, with the added advantage of not scuffing enamel.

IPR Tips

IPR™ Tips allow interproximal reduction safely and with optimal control.

The tip thickness sets the prescribed amount of reduction. They are double-sided, except for the 0.25 mm tips which come with abrasive on the left or right side. Like all Profin tips, they do not harm soft (gingival) tissue.



EVA Tips



EVA™ Tips are made of plastic. They are used for cleaning and polishing, implants, prosthetics, restorations and more, in interproximal and other hard-to-reach spaces. They come in two types for use with polishing paste.


Lamineer Tip Kits with Stand

10 Lamineer Tips w/ Stand

Item No: LTA

Includes:

- 10 Safe-Sided Lamineer Tips
 - Diamond Short Spatulas: LTA-S100, LTA-S50, LTA-S15
 - Diamond Standard Spatulas: LTA-150, LTA-75, LTA-30
 - Diamond Shoulder/Bevel Tip: LTA-C50
 - Diamond Concave Tip: LTA-E50
 - Wolfram Tungsten Spatulas: LTA-S36, LTA-39
- 1 Stand w/ Color Coded Chart



10 Lamineer Tips w/ Short Spatulas

Item No: LTA-S

Includes:

- 10 Safe-Sided Lamineer Tips w/ Short Spatulas
 - Diamond Short Spatulas: LTA-S150, LTA-S100, LTA-S75, LTA-S50, LTA-S30, LTA-S15
 - Diamond Shoulder/Bevel Tip: LTA-C50
 - Diamond Concave Tip: LTA-E50
 - Wolfram Tungsten Spatulas: LTA-S36, LTA-39
- 1 Stand w/ Color Coded Chart



20 Lamineer Tips w/ Stand

Item No: LTA-U

Includes:

- 20 Safe-Sided Lamineer Tips
 - Diamond Short Spatulas: LTA-S150, LTA-S100, LTA-S75, LTA-S50 (2pcs), LTA-30, LTA-S15
 - Diamond Standard Spatulas: LTA-150, LTA-100, LTA-75, LTA-50 (2pcs), LTA-30, LTA-15
 - Diamond Shoulder/Bevel Tip: LTA-C50
 - Diamond Concave Tip: LTA-E50
 - Wolfram Tungsten Spatulas: LTA-S36 (2pcs), LTA-39 (2pcs)
- 1 Stand w/ Color Coded Chart



Lamineer Tip Refills

Diamond Standard Spatula Tips

Diamond Standard Spatulas are for reshaping and gross reduction of overhangs of all restorative materials and enamel.



Diamond Spatula Tip Refills

	Item No.	Description
	LTA-150/2	2 Spatulas, 150 Micron, Safe-Sided
	LTA-100/2	2 Spatulas, 100 Micron, Safe-Sided
	LTA-75/2	2 Spatulas, 75 Micron, Safe-Sided
	LTA-50/2	2 Spatulas, 50 Micron, Safe-Sided
	LTA-30/2	2 Spatulas, 30 Micron, Safe-Sided
	LTA-15/2	2 Spatulas, 15 Micron, Safe-Sided

Wolfram Tungsten Standard Spatula Tips

Wolfram Tungsten Standard Spatulas are for finishing and polishing of all restorative materials and enamel. They do not scuff healthy enamel.



Wolfram Tungsten Spatula Tip Refills

	Item No.	Description
	LTA-39/2	2 Spatulas, Safe-Sided, Ultra-thin
	LTA-38/2	2 Spatulas, Two-Sided, Ultra-thin

Lamineer Tip Refills

Diamond Short Knife-Edge Spatula Tips

Diamond short knife-edge spatulas are for reshaping and gross reduction of overhangs of all restorative materials and enamel.

They are the same as the standard spatula but shorter and narrower with knife edge for better interproximal access.



Diamond Short Knife-Edge Spatula Tip Refills

	Item No.	Description
	LTA-S150/2	2 Short Knife-Edge Spatulas, 150 Micron, Safe-Sided
	LTA-S100/2	2 Short Knife-Edge Spatulas, 100 Micron, Safe-Sided
	LTA-S75/2	2 Short Knife-Edge Spatulas, 75 Micron, Safe-Sided
	LTA-S50/2	2 Short Knife-Edge Spatulas, 50 Micron, Safe-Sided
	LTA-S30/2	2 Short Knife-Edge Spatulas, 30 Micron, Safe-Sided
	LTA-S15/2	2 Short Knife-Edge Spatulas, 15 Micron, Safe-Sided

Wolfram Tungsten Short Knife-Edge Spatula Tips

Wolfram Tungsten short knife-edged spatulas are for finishing and polishing of all restorative materials and enamel. They do not scuff healthy enamel. They are the same as standard spatula but shorter and narrower with knife edge for better interproximal access.



Wolfram Tungsten Short Knife-Edge Tip Refills

	Item No.	Description
	LTA-S36/2	2 Short Knife-Edge Spatulas, Safe-Sided

Lamineer Tip Refills

Diamond Concave Tips

Diamond Concave Tips are for fine protrusive and lateral adjustments on concave surfaces, e.g. cingulum and lingual fossa areas.




Diamond Concave Tip Refills

	Item No.	Description
	LTA-E50/2	2 Diamond Concave Tips, 50 Micron

Diamond Shoulder/Bevel Tips

Diamond Shoulder/Bevel Tip is for shaping precise inclines and angles of shoulder and bevels.




Diamond Shoulder/Bevel Tip Refills

	Item No.	Description
	LTA-C50/2	2 Diamond Shoulder/Bevel Tips, 50 Micron

Diamond Safe Sided Tips

Diamond safe sided tips are for stripping and shapping contacting surgaces.



Diamond Safe Sided Tip Refills

	Item No.	Description
	LTA-L75/2	2 Diamond Left Abrasive Side Tips, 75 Micron
	LTA-L50/2	2 Diamond Left Abrasive Side Tips, 50 Micron
	LTA-R75/2	2 Diamond Right Abrasive Side Tips, 75 Micron
	LTA-R50/2	2 Diamond Right Abrasive Side Tips, 50 Micron

Lamineer Tip Refills

Diamond Paddle Tips

Diamond Paddle Tips are for refinishing prepared surfaces of all restorative materials and enamel. They remove grooves and undercuts.



Diamond Paddle Tip Refills

	Item No.	Description
	LTA-D150/2	2 Diamond Paddle Tips, 150 Microns
	LTA-D100/2	2 Diamond Paddle Tips, 100 Microns
	LTA-D50/2	2 Diamond Paddle Tips, 50 Microns
	LTA-D30/2	2 Diamond Paddle Tips, 30 Microns

Diamond Saw Blade Tips

Diamond Saw Blade Tips are safe sided and used for opening incisal and gingival embrasures.




Diamond Saw Blade Tip Refills








	Item No.	Description
	LTA-T50/2	2 Diamond Saw Blade Tips, 50 Microns
	LTA-T30/2	2 Diamond Saw Blade Tips, 30 Microns

IPR Tip Refills

IPR Tip color coding corresponds to the thickness of the blade.
They are abrasive on both sides other than the P25L and P25R
which are abrasive on only one side.



IPR Tip Refills



	Item No.	Description
	LTA-P5/2	2 IPR Tips, 0.5mm thick, 50 microns
	LTA-P4/2	2 IPR Tips, 0.4mm thick, 50 microns
	LTA-P3/2	2 IPR Tips, 0.3mm thick, 50 microns
	LTA-P25L/2	2 IPR Tips, 0.25mm thick, abrasive left side, 50 microns
	LTA-P25R/2	2 IPR Tips, 0.25mm thick, abrasive right side, 50 microns
	LTA-P2/2	2 IPR Tips, 0.2mm thick, 40 microns
	LTA-P1/2	2 IPR Tips, 0.1mm thick, 15 microns

EVA Tip Refills

EVA tips are made of plastic. and used for cleaning and
polishing with polishing paste.



EVA Tip Refills

	Item No.	Description
	EVA-2000/100	100 Tips, V-Shaped with collapsible sides
	EVA-5000/100	100 Tips, Paddle with round tip